

Pi-Clear nucleic acid purification kit - DNA/RNA viral

ATTENTION: Before starting the assay, read the item "Special Care" in "User Guide".

Spin Column – PB1	50 units (spin column + collection tube)
Proteinase K - PK	300 µL
Carrier - CR	350 µg (freeze dried reagent. Resuspend in 350 µL of PB3)
Lysis Buffer – TL	15 mL
Wash Buffer 1 – TP1	50 mL
Wash Buffer 2 – TP2	OBS: 48 mL of Ethanol 96-100 % must be added Final Volume: 60 mL
Recovery tube (1,5 mL) – PB2	50 units
Nuclease-free Water – PB3	15 mL

Storage: PB1, TL, TP1, TP2, PB2 and PB3 must be stored at room temperature (15 - 25 °C); PK and CR (after resuspension in PB3) must be stored at -20 °C.


Required materials: Ethanol 96-100% and extra DNase/RNase-free microtubes.

Validated procedure: sputum, liquor, swab, semen, whole blood, blood plasma and blood serum.

Procedure

A) Sample Lysis


- 1 Add 200 µL of sample into a clean, sterile and DNase/RNase free microtube;
Obs: prepare sample as described in "Instruction for Use";
- 2 Add 5 µL Proteinase K (PK) and vortex;
- 3 Add 200 µL Lysis Buffer (TL) and vortex.



Nucleic acid Extraction and Purification


B) Binding

- 4** Add 5,6 μL Carrier (CR) and vortex;
- 5** Incubate at room temperature for 3 minutes;
- 6** Add 200 μL Ethanol (96-100 %) and vortex;
- 7** Incubate at room temperature for 5 minutes.
- 8** Transfer up total volume of lysed sample to spin column with the collection tube (PB1);
- 9** Centrifuge at 4.000 xg for 3 minutes;
- 10** Discard the flow-through.




C) Washing

- 11** Add 700 μL Wash Buffer 1 (TP1);
- 12** Centrifuge at 11.000 x g for 1 minute;
- 13** Discard the flow-through;
- 14** Add 500 μL Wash Buffer (TP2);
- 15** Centrifuge at 11.000 x g for 1 minute;
- 16** Discard the flow-through;
- 17** Repeat steps 14-16 once;
- 18** Centrifuge at 12.000 x g for 2 minutes to dry the membrane with bound DNA/RNA;
- 19** Discard the collection tube and insert the spin column into a recovery tube (PB2).



D) Elution

- 20** Add 30-100 μL Nuclease-free Water (PB3) to the center of spin column;
- 21** Incubate at room temperature for 3 minutes;
- 22** Centrifuge at 15.000 x g for 3 minutes and discard the spin column;
- 23** Store your purified DNA/RNA at -20 °C. For long-term storage, keep the purified DNA/RNA at -70 °C.



Access the complete "User Guide"